

# 118 Helpsheet

---



*World Leaders in Computer Controlled Testing  
Systems for Geotechnical Engineers and Geologists*

Hardware

IEEE Interface

Selecting interface card for use with GDS hardware

---

## 1. Introduction

---

GDSLAB Windows software uses the National Instruments PCI-GPIB interface card to communicate with GDS GPIB hardware using NI-488.2 communication. The interface card can be purchased from GDS or directly from your local National Instruments supplier.

## 2 N.I. GPIB Interface Card Part Number

---

To order your interface card directly from your local National Instruments supplier quote the following part number:

778032-01 PCI-GPIB, NI-488.2M for Windows 2000/XP/ME/98

777158-01 PCI-GPIB and NI-488.2 for win 95

This interface card takes up one of the PCI slots inside your computer.